

Brighton Pediatric Early Warning Score

	0	1	2	3
Behavior	Playing Appropriate	Sleeping	Irritable	Lethargic &/OR Confused &/OR Reduced response to pain
Respiratory	Within normal parameters, no recession OR tracheal tug	10 above normal parameters, using accessory muscles &/OR 30+ % FiO2 or 4+ liters/min	>20 above normal parameters, recessing/retractions, tracheal tug &/OR 40+ % FiO2 or 6+ liters/min	5 below normal parameters with sterna recessions/retractions, tracheal tug or grunting &/OR 50% FiO2 or 8+ liters/min
Cardiovascular	Pink &/OR Capillary refill 1–2 seconds	Pale &/OR Capillary refill 3 seconds	Grey &/OR Capillary refill 4 seconds, Tachycardia of 20 above normal rate	Grey & mottled or capillary refill 5 seconds or above OR Tachycardia of 30 above normal rate of bradycardia
Bronchodilators every 15 minutes &/OR Persistent vomiting following surgery (2 points each)				
TOTAL PEWS SCORE:				

Client: CHBC

Date: 2 November 2015 1:15 PM

File Name: 21012223-P_CHBC-Lanyard-Brighton_SBAR

Target: Lanyards

Actual Size: 95mm x 50mm

Revisions: 0

Colours: 4c

Operator: CR

SBAR Communication Tool

S **Situation:** This is _____, and I am the nurse for patient _____ on pediatrics. I am concerned about the patient because _____ (give concise description — Pain, VS change, status change, increasing PEWS score).

B **Background:** Patient was admitted on _____ with _____. Recent procedures include (OR, tests, investigations) _____. The child's condition has changed in the last _____ minutes. The latest vital signs are HR, RR, BP, O2 sats, Temp, Pain Score, PEWS

A **Assessment:** I think the problem is _____. I have given (O2, analgesia, stopped infusion) _____. The patient seems to be unstable, deteriorating, needs intervention.

R **Recommendation:** I need you to:

- Come and see the patient now
- Call the MRP or pediatrician
- Order tests or labs (CXR, Labs)

Is there anything I need to do now (monitoring, fluids,)
End call by confirming plan of care



Client:	CHBC	Date:	2 November 2015 1:15 PM
File Name:	21012223-P_CHBC-Lanyard-Brighton_SBAR	Target:	Lanyards
Actual Size:	95mm x 50mm	Revisions:	0
Colours:	4c	Operator:	CR